
VP-D1

CE 0482

innoGIO

GIÒvital

COMPRESSOR NEBULIZER

USER'S MANUAL

Manual Version: 1.1.
Date of Issue: AUG 2019

Contents

Before Using the Unit

Notes on Safety	3
Product Introduction	6
Product Contents	7
Name of the Parts	8

How to Assemble

How to Prepare the Power Source	9
How to Fill the Medication	9
Assemble the Medicine Cup	10
Assemble the Air Tube	10

How to Use the Unit

How to Use the Unit	11
---------------------------	----

How to Clean, Disinfect and maintain the device

How to Clean the Unit after Inhalation	14
How to Disinfect	14
How to Clean the Air Filter	15
Maintenance of the Device	16

Troubleshooting

Troubleshooting	18
Specifications	20
Guidance and Manufacturer's Declaration	21

Notes on Safety

Be sure to read the notes on safety before using the unit.

The warning signs and the sample icons shown here are listed for you use this product safely and correctly as well as to prevent the risk and the damage to you and others from happening. The icons and meanings are as follows.

Warning

- For the type, dose, and regimen of the medication, be sure to follow the instructions of a doctor.
- When you use the unit for the first time after purchasing it or after not using it for a long period of time, be sure to clean and disinfect Medicine Cup, Inhalation Mask Mouthpiece and Nosepiece (optional).
- Clean and disinfect Medicine Cup, Inhalation Mask, and Mouthpiece after each use.
- Be sure that the patient is an intended operator.
- Don't maintain the device while it is in use with the patient.
- Clean and disinfect the Medicine Cup, Mouthpiece, Inhalation Mask and Nosepiece (optional) before using them for the first time after purchase.
- Be sure to dry the cleaned and disinfected parts promptly, and store them in a clean place.
- Do not place or attempt to dry the device, components or any of the nebulizer parts in a microwave oven.
- Keep the device out of the reach of unsupervised infants and children. The device may contain small parts that can be swallowed, moreover strangulation due to cables and hoses may occur.
- Don't use accessories such as medicine cup, mask, nosepiece, mouthpiece what is not provided by the supplier described in the brochure.
- Do not wash the Main Unit and Power Plug with water or hot water. Do not immerse the Main Unit in water or other liquid.
- Do not cover the compressor with a blanket, towel, or any other type of cover during use. This could result in the compressor overheating or malfunctioning.
- Do not use the device where the device may be exposed to flammable gas or vapors. Pentamidine is not an approved medication for use with this device.
- Always dispose of any remaining medication in the medicine cup after each use.

- Do not leave the device or its parts where it will be exposed to extreme temperatures or changes in humidity, such as leaving the device in a vehicle during warm or hot months or where it will be exposed to direct sunlight.
- Do not use or store the device where it may be exposed to noxious fumes or volatile substances. Make sure that the Medicine Cup is clean before use.
- Do not use in anaesthetic or ventilator breathing circuits.
- Not suitable for use in the presence of flammable anesthetic mixture with air or oxygen or nitrous oxide.
- Risk of Electrical Shock:
Do not use the compressor (main unit) or Power Plug while they are wet. Do not plug or unplug the Power Plug into the electrical outlet with wet hands. Do not use or store the device in humid locations, such as a bathroom. Use the device within the operating temperature and humidity. Do not operate the device with a damaged power cord or plug.
- Maintenance and Storage
Do not leave the cleaning solution in the nebulizer parts. Rinse the nebulizer parts with clean hot tap water after disinfecting.
- Store the device and the components in a clean, safe location.
- If you have any problems with this device, such as setting up, maintaining or using, please contact with service personnel of INNOGIO. Don't open or repair the device by yourself.

Caution

- Do not inhale water in the Medicine Cup.
- Do not spill water or other liquids on the compressor and Power Plug. These parts are not waterproof. If liquid spills on these parts, immediately unplug the Power Plug and wipe off the liquid with gauze or other soft absorbent material.
- Do not drop or apply strong shock to the Main Unit or Medicine Cup.
- Do not disassemble, repair, or modify the unit.
- Provide close supervision when this device is used by, on, or near infants, children or compromised individuals.
- If the device is used continuously, the service life of the device may be shortened.
- Limit use to 25 minutes at a time, and allow a 40 minutes interval before using the device again.
- Do not insert any object into the compressor.

- Make sure the Air Filter is clean. If the Air Filter has changed color or has been used for more than 60 days, replace it with a new one.
- Make sure the Medicine Cup is correctly assembled, the Air Filter is properly installed, and the Air Tube is correctly connected to the compressor and the Medicine Cup. Air may leak from the Air Tube during use if not securely connected.
- Do not use the device if the Air Tube is bent.
- Do not add more than 8 ml of medication to the Medicine Cup.
- Do not operate the device at temperatures greater than +40°C (+104°F).
- Do not tilt the Medicine Cup so the angle of the kit is greater than 45°.
- Medication may flow into the mouth. Do not shake the Medicine Cup while using the device.
- Use only InnoGIO authorized parts and accessories. Parts and accessories not approved for use with the device do not perform the expected specification or it may damage the unit.
- Changes or modification not approved by InnoGIO will void the user warranty.
- To avoid injury to the nose mucosa, do not squeeze the optional Nosepiece into the back of the nose.
- When using this device, there will be some noise and vibration caused by the pump in the compressor.
- There will also be some noise caused by the emission of compressed air from the Medicine Cup.
- This is normal and does not indicate a malfunction.
- When sterilizing parts by boiling, make sure that the container does not boil dry.
- Do not use the device while sleeping or if drowsy.
- Remove the Power Plug from the device after use.
- Unplug the Power Plug from the electrical outlet before cleaning the device.
- Do not leave the Medicine Cup with medication after use.
- Do not store the Air Tube with moisture or medication remaining in the Air Tube. This could result in infection as a result of bacteria.
- Do not wrap the Power Plug around the compressor (main unit).
- People who is not trained with specified skills are not permitted to maintain the device.
- If you have any problems with this device, such as setting up, maintaining or using, please contact with service personnel of InnoGIO. Don't open or repair the device by yourself.

Product Introduction

This product is intended for the treatment of respiratory diseases such as asthma allergies etc. Working principle: The compressor forces air to the nebulizer. When the air enters the nebulizer, it converts the prescribed medication into aerosol of microscopic droplets that can easily be inhaled.

The meaning of the symbols used in this instruction brochure.

Symbols

Symbol for caution. This is a safety symbol that is used to highlight that there are specific warnings or precautions associated with the devices that are not otherwise found on the label. The symbol can also be used to mean "attention, see Instructions for Use."

Manufacturer

The symbol indicates that the device should be sent to the special agencies according to local regulations for separate collection after its useful life.

Symbol for "THE OPERATION GUIDE MUST BE READ"; Indicates the need for the user to consult the instructions for use.

This symbol indicates that the device includes IEC 60601-1 Type BF Applied Part.

IPX1

Protection against ingress of water, it means the device could protected against vertical falling water drop 15°.

Indicates the date when the medical device was manufactured.

SN

This symbol shall be accompanied by the manufacturer's serial number.

LOT

Symbol for lot number; Indicates the manufacturer's batch code so that the batch or lot can be identified.

EC REP

Indicates the Authorized Community.


Symbol for CE Mark. This symbol certifies that a product has met European Union consumer safety, health, or environmental requirements.

Product contents

Please check the product.

The following items are contained in the box. If any item is missing, please contact the store where you purchased the product or the nearest InnoGIO dealer.

1. Compressor
(Main Unit)

2. Medicine Cup

3. Inhalation Mask
(Adult)

4. Inhalation Mask
(Children)

5. Mouthpiece

6. Nosepiece
(Optional)

7. Air Tube

8. Air Filters
(Optional)


9. DC Adaptor

To start inhalation, you must assemble the unit. Please read the Page 9-11 and understand the points for assembling before starting to use the unit.


(Refer to Pages 9-11)

Name of the Parts

Name of each part

- | | |
|---------------------------|------------------------------|
| 1. Compressor (Main Unit) | 8. Air Filter Cover |
| 2. Air outlet | 9. Footpad |
| 3. Air Tube | 10. DC adaptor |
| 4. Power Switch | 11. Medicine cup cover |
| 5. DC jack | 12. Medicine cup upper shell |
| 6. Medicine cup | 13. Nebulizer slice |
| 7. Medicine cup Plug | 14. Medicine cup |

Name of accessories

- a) Medicine Cup
- b) Inhalation Mask (child/adult)
- c) Mouthpiece
- d) Nosepiece (optional)
- e) Air Tube
- f) Air Filters (optional)
- g) DC Adaptor

How to Assemble


Make sure of that the mouthpiece, nosepiece, mask and air filter be clean before use.

- Clean and disinfect the Medicine Cup, mouthpiece, nosepiece or optional masks before using them for the first time after purchase (refer to page 14).
- If the Unit has not been used for a long period of time, or if more than one person uses the same device, clean and disinfect the Medicine Cup, mouthpiece, nosepiece or optional masks before using them (refer to page 14).
- Make sure the air filter is clean. If the air filter has changed color or has been used on average for more than 60 days, replace it with a new one (refer to page 15).

1. Prepare the Power Source

Connect the DC adaptor to the DC Jack of the Main Unit, and then to the power outlet as shown in the Figure.

Don't operate with wet hand!


2. Fill in the Medication

Lift up the Mouthpiece, or Nosepiece or Mask from the Medicine Cup, and unscrew the Nebulizer Upper Shell.

Caution: Be sure not to knock the baler off to the ground.


The baler is fixed on the Nebulizer Upper Shell by two hooks, and can be pulled out with some power.

Fill in the correct amount of prescribed medication in the Nebulizer Bottom Shell.


3. Assemble to the Medicine Cup

A. Screw the Nebulizer Upper Shell to the Nebulizer Bottom Shell clockwise until it is securely closed.


B. Assemble the desired inhalation accessories.

Assemble the Mouthpiece

Assemble the Nosepiece


Assemble the Mask

4. Assemble the Air Tube


Step 1: Twist the Air Tube plug and push it firmly into the Air Tube Connector of the compressor.

Step 2: Twist the Air Tube plug slightly and push it firmly into the Air Tube connector Air Tube connector on the bottom of the Medicine Cup.

Caution: Keep the Medicine Cup upright when assembling the Air Tube. Otherwise medication may spill out the Medicine Cup.
Caution: Don't use the Unit if the Air Tube is bent.

Use the Medicine Cup Holder as a temporary holder for the Medicine Cup.


How to Use the Unit

1. Hold the Medicine Cup as shown in the Figure


Correct Posture

Caution: Do not tilt the Medicine Cup greater than 45°. Otherwise medication may flow into the mouth.

2. Turn on the Unit

Turn the Power Switch to the on (I) position. As the compressor starts, nebulization begins. (To stop nebulization, turn the power switch to the off (O) position.)

3. Inhale medication as instructed by physician

Difference among the inhalation accessories as following.

- **Mouthpiece**

Insert the Mouthpiece into your mouth.
Inhale medication normally.
Exhale normally through the Mouthpiece.

- **Mask**

Place the Mask over the nose and mouth.
Pull the elastic strap over the head.
Gently pull on the strap to secure the mask
over the nose and mouth. Inhale the medication.
Exhale normally through the mask.


- **Nosepiece (Optional)**

Inhale the medication through the nosepiece
as illustrated. Exhale through the mouth.

4. Turn off the Unit

Turn the Power Switch
to the off (O) position.
As the compressor stops,
nebulization stops.

5. Disconnect the Air Tube from the Medicine Cup

- Hold the Air Tube Plug and gently pull it down.
- Check the Air tube. No condensation or moisture should remain in the Air Tube.

If any condensation or moisture remains in the Air tube, remove the moisture from the air tube. Follow the directions below:


1. Make sure the Air Tube is still connected to the air tube connector of the compressor.
2. Turn on the Unit, and the compressor will start and pump air through the air tube to expel the moisture.
3. Turn off the Unit.

6. Disconnect the air tube from the compressor

- Hold the Air Tube Plug and gently pull the Air Plug off the Air Tube Connector on the upper side of the compressor.

7. Unplug the DC Adapter from the power outlet

- Unplug the DC Adaptor from the power outlet and DC Jack the main unit.


How to Clean the Unit after Inhalation

Be sure to clean and store the unit after inhalation.

1. Remove the inhalation accessory from the Medicine Cup.
2. Disconnect the Air Tube from the Medicine Cup.
3. Gently screw off the Nebulizer Upper Shell.
4. Discard the remaining medication.
5. Wash the parts suitably in water.
6. Hand dry or air dry in a clean environment using a soft, clean, lint-free cloth.
7. Assemble the Medicine Cup and store it in a dry bag.

Caution: Clean the casing of the main unit by using a soft cloth moistened with water or mild detergent. Do not use abrasive cleaners.


How to Disinfect

Disinfect the Medicine Cup and Mouthpiece or nosepiece after the last treatment of the day.

- The method as described below:

Use a commercially available disinfectant with the instructions provided by the disinfectant manufacturer, it is applicable to the masks' disinfection.



- Submerge the parts in the disinfectant solution for the specified period.
- Remove the parts and discard the solution.
- Rinse the parts with clean hot tap water, shake off excess water and allow to air dry in a clean environment.

- Never dry the parts in the microwave oven.
- Never use an autoclave, EOG gas disinfection or low temperature plasma sterilizer to disinfect the device.

How to clear the Air Filter

1. Pull the Air Filter cover to remove from the back side of the compressor.
2. Remove the dirty Air Filter.
CAUTION
Do not wash or clean the Air Filter. It may cause blockage.
Do not use cotton or other material to substitute the Air Filter.
3. Insert a new Air Filter.
CAUTION
Before inserting the new Air Filter make sure the Air Filter is clean.
Do not operate the device without the Air Filter.
Use only the InnoGIO Air Filter.
4. Put the Air Filter cover back to the compressor.

NOTE:

If the Air Filter has changed color or has been used for more than 60 days, replace it with a new one. If water or medication has spilled on the Air Filter, do not wash or clean the Air Filter. It may cause blockage. Do not use cotton or other material to substitute the Air Filter. Before inserting the new Air Filter make sure the Air Filter is clean. Do not operate the device without the Air Filter. Use only the InnoGIO Air Filter.

Maintenance of the Device

USE

This unit is a medical instrument. Be sure to follow the instructions of a doctor and use the unit correctly.

The nebulization characteristics of this unit differ by the properties of medication. Especially with the use of a medication having high surface activity or viscosity such as medication solubilizing agent or expectorant, the nebulization rate may be reduced. The nebulization rate may also be reduced when the temperature of medication is low.

Do not use the device if it is damaged in any way. The continuous use of a damaged unit may cause injury, improper results, or serious danger.

When not in use, store the device in a dry room and protect it against extreme moisture, heat, lint, dust and direct sunlight. Never place any heavy objects on the storage case.

Provide close supervision when the device is used by, on, or near infants, children or compromised individuals. No modification of this equipment is allowed.

Inhalation Mask is recommended when the main unit is used by children less than 5 years old. Keep the device out of the reach of unsupervised infants and children. The mouthpiece or other small parts may be swallowed.

Metallic sound may be generated, but it is not a trouble.

The device shall be disconnected from the power source after use.

The mouthpiece and inhalation mask contact user directly, please keep it clean.

STORAGE

DO NOT store the device in locations exposed to high temperatures (over 50°C), low temperatures (below -20°C), high relative humidity (over 85%) or excessive amounts of dust. Never place any heavy objects on the storage case.

Store the device and component in a clean, safe location. Keep the unit out of reach of small children.

MAINTENANCE

DO NOT drop the device or subject it to other shocks or vibration.

Never dry the parts in the microwave oven.

Be noted that parts other than mouthpiece and inhalation mask cannot be boiled!

Do not carry or leave the nebulizer with medication in the medicine cup.

*If you have any problems with this device, such as setting up, maintaining or using, please contact with service personnel to InnoGIO.
Don't open or repair the device by yourself.*

Troubleshooting problems

Trouble	Possible Cause	How to Correct
No power to the device when the Power Switch is on.	The DC adapter is not plugged into the main unit.	Turn the Power Switch off. Plug the DC adapter into the DC jack. Turn the device on.
No nebulization or low nebulization rate when the power is on.	No medication in the Medication Cup. Too much or too little medication in the Medication Cup.	Add the correct amount of prescribed medication to the Medication Cup.
	The Baffler is not attached to the Medicine Cup Upper Shell or incorrectly positioned.	Make sure the Baffler is correctly attached to the Medicine Cup Upper Shell.
	The Medicine Cup is not correctly assembled.	Make sure the Medicine Cup is correctly assembled and the inhalation accessory is correctly attached.
	The nozzle is blocked.	Clean and disinfect the Medicine Cup to remove the blockage
	The Medicine Cup is tilted at an incorrect angle.	Hold the Medicine Cup correctly. Do not tilt the Medicine Cup so the angle of the kit is greater than 45 degrees.
	The Air Tube is incorrectly attached.	Make sure the Air Tube is correctly attached to the compressor and the Medicine Cup.
	The Air Tube is folded or damaged. The Air Tube is blocked.	Make sure the Air Tube is not folded, kinked or bent. Inspect the Air Tube for any damage. Replace the Air Tube if damaged
	The Air Filter is dirty	Replace the Air Filter with a new clean Air Filter

Troubleshooting problems

The device is very hot.	The compressor is covered.	Do not cover the compressor with any type of cover during use.
	Operating continuously over 25 minutes.	Limit use to 25 minutes at a time and allow a 40 minutes interval before using the device again.
The device is abnormally loud.	The Air Filter cover is incorrectly attached.	Attach the Air Filter cover correctly. Make sure the Air Filter cover is not blocked.

Note: If the unit does not nebulize normally after taking the above-mentioned procedure, contact the store where you purchased the unit or the nearest InnoGIO dealer.

Specifications

Product name	InnoGIO GIO Vital Compressor Nebulizer
Model	VP-D1
Applied Part	Mask, Mouthpiece, Nosepiece (optional), Type BF
Power source	AC 100V ~ 240V,50Hz/60Hz
Working Environment	Temperature: 10 to 40 °C Humidity: 30 to 85% RH Atmospheric Pressure: 700 – 1060hPa
Storage Environment	Temperature: -20 to 50 °C Humidity: 30 to 85% RH Atmospheric Pressure: 500 – 1060hPa
Classification of installation and use	Main unit: portable, Applied part: Hand-held Non-AP, APG portable
Device type	Equipment
Power consumption	adaptor: I/P:100-240VAC, 50/60Hz, 0.35MAX. O/P: 12VDC, 1A. main unit: I/P: 12VDC, 1A
Sound Level	≡52dB
Capacity of Medicine Cup	≡8 ml
Nebulization Rate	±0.20 ml / min
Particle Size	MMAD ≡ 5 μm
Dimensions	Model VP-D1 123mm(L)*81mm(W)*53.9mm(H)
Weight (not including battery)	Model VP-D1 Approximately 228g (not including Medicine Cup and Air Tube)
Pollution	Degree 2
Voltage Degree	CATEGORY I
Altitude Category	≡3000 m
Protection against harmful ingress of water or particulate matter	IPX1
Free Flow Rate	4L/min–7L/min
Nebulization Pressure	20Kpa–80Kpa (2.9psi-11.7psi)

Nebulization rate is lowered at low temperature.

Specifications and appearance may be changed without prior notice for improvement.

Guidance and Manufacturer's Declaration

1. This equipment needs to be installed and put into service in accordance with the information provided in the ACCOMPANYING DOCUMENTS;
2. IMMUNITY TEST LEVELS for basic safety and essential performance of ME equipment and ME systems should be chosen based on a high probability of maintaining basic safety and essential performance, and shall be according to the professional healthcare facility environment, home healthcare environment, and special environment, based on the locations of intended use.
3. HOME HEALTHCARE ENVIRONMENT is dwelling place in which a patient lives or other places where patients are present, excluding professional healthcare facility environments where operators with medical training are continually available when patients are present. Such as schools, outdoors, domiciles, vehicles hotels and pensions.
4. Please dispose the parts of the device once it is out of use according to local authorities.
5. MANUFACTURER will provide circuit diagrams, component part lists, descriptions, calibration instructions to assist to SERVICE PERSONNEL in parts repair.

EXAMPLE: As indicated in Table 6 of IEC 60601-1-2:2014 for ME EQUIPMENT, a typical cell phone with a maximum output power of 2 W yields $d = 3.3$ m at an IMMUNITY LEVEL of 3 V/m.

A1 Electromagnetic Emissions – For all EQUIPMENT and SYSTEMS

Guidance and manufacture's declaration-electromagnetic emission

The InnoGIO Compressor Nebulizer, model GIOVital VP-D1 is intended for use in the electromagnetic environment specified below; The customer or the user of the InnoGIO Compressor Nebulizer, model VP-D1 should assure that it is used in such and environment.

Emission test	Compliance	Electromagnetic environment-guidance
Rf emissions CISPR 11	Group 1	The InnoGIO Compressor Nebulizer, model VP- D1 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Rf emissions CISPR 11	Class B	The InnoGIO Compressor Nebulizer, model VP-D1 is suitable for use in all establishments, other than domestic and those directly connected to the public low-voltage power supply network that supplies building used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

A2 Electromagnetic Immunity – For Home Healthcare Environment EQUIPMENT and SYSTEMS

Guidance and manufacture's declaration – electromagnetic immunity

The InnoGIO Compressor Nebulizer, model GIOVital VP-01 is intended for use in the electromagnetic environment specified below. The customer or the user of InnoGIO Compressor Nebulizer, model VP-D1, should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment -guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	±8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.

Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±2KV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line to line ±2 kV line to ground	±1 kV line to line ±2 kV line to ground	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30A/m	30A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the InnoGIO Compressor Nebulizer, model VP-D1 requires continued operation during power mains interruptions, it is recommended that the VAPO Compressor Nebulizer, model VP-D1 be powered from an uninterruptible power supply.
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHZ 10 V/m 80 MHZ to 2.7 GHz	3 Vrms 10 V/m	The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation: $E=6/d \sqrt{P}$ Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.

NOTE UT is the a.c. mains voltage prior to application of the test level.

Importer:

Innogio Sp. z o.o., Ostródzka street 74H, 03-289 Warsaw, Poland

innogio


Authorized European Representative:

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MedNet GmbH